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OPERATION

1. Attach the amplifier to the handset tightly as shown.

2. Turn the amplifier POWER INDICATOR ON. The power indicator will light up. Adjust

the volume to suit your hearing.

- 3. If you are using this amplifier together with a hearing aid equipped with a telecoil, switch the hearing aid to "T" or "Loop" position. The amplifier functions as an inductive coupler with the telecoil to minimize the background noises and thus enables you to hear more clearly.
- 4. If you experience "feedback" (a' high-pitched squeal), you

- may have the volume up too high for your telephone. Correct this by reducing the volume.
- 5. If the voice sounds low or distorted, or if the power indicator is dim or fails to light up, your batteries may be low. Check or replace the batteries.
- 6. When finished, turn dial to "OFF" and remove the amplifier from the handset.

WARRANTY

Ameriphone will repair your Portable Telephone Amplifier for six months from the date of purchase if the unit is defective in workmanship or materials. (Batteries are not included in this warranty). Proof of purchase is required.

All applicable implied warranties, including the implied warranty of merchantability and fitness for a particular purpose, given to you by law are hereby limited in

duration of this warranty. **UNDER NO CIRCUMSTÁNCES** WILL AMERIPHONE BE LIABLE FOR ANY INCIDENTAL OR **CONSEQUENTIAL DAMAGES.**

You may also have other warranty rights which vary from state to state.

After the six month warranty period, Ameriphone will make repairs at a minimal charge. Contact our service department for details

QUALITY PRODUCTS

Ameriphone offers an extensive line of quality communications products. Call for a FREE catalogue, today.

ameriphone

12082 Western Ave., Garden Grove, CA 92841 (714) 897-0808 VOICE (714) 897-1111 TTY (714) 897-4703 FAX

© AMERIPHONE, Inc. 5210-0261 🛕 5/00 M02798 ENGLISH

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9/OH9199/AG